

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD94

Numéro de catalogue: **CL750-65240**



Informations de base

Numéro de catalogue:

CL750-65240

Taille:

100ug, 0.5 mg/ml

Hôte:

Rat

Isotype:

IgG2a, kappa

Numéro d'acquisition GenBank:

BC120854

Identification du gène (NCBI):

16643

Nom complet:

killer cell lectin-like receptor,
subfamily D, member 1

Méthode de purification:

Purification par affinité

CloneNo.:

18D3

Excitation/Emission maxima
wavelengths:

755 nm / 780 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

souris

Informations générales

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ mouse splenocytes were surface stained with FITC Anti-Mouse CD49b (HMA2) (FITC-65245, Clone: HMA2) and 0.5 ug CoraLite® Plus 750 Anti-Mouse CD94 (CL750-65240, Clone: 18D3) or CoraLite® Plus 750 Rat IgG2a isotype control. Cells were not fixed.