

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD49b

Numéro de catalogue: **CL750-65252**



Informations de base

Numéro de catalogue:

CL750-65252

Taille:

100 tests, 5 µl/test

Hôte:

Mouse

Isotype:

IgG1, kappa

Numéro d'acquisition GenBank:

NM_002203

Identification du gène (NCBI):

3673

Nom complet:

integrin, alpha 2 (CD49B, alpha 2 subunit of VLA-2 receptor)

MW calculé

126 kDa

Méthode de purification:

Purification par affinité

CloneNo.:

AK7

Excitation/Emission maxima wavelengths:

755 nm / 780 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

Humain

Informations générales

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ human PBMCs were surface stained with 5 ul CoraLite® Plus 750 Anti-Human CD49b (CL750-65252, Clone: AK7) or unstained. Cells were not fixed. Monocytes were gated.