

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD49b

Numéro de catalogue: FITC-65245



Informations de base

Numéro de catalogue:

FITC-65245

Taille:

100ug , 500 µg/ml

Hôte:

Armenian Hamster

Isotype:

IgG

Numéro d'acquisition GenBank:

BC065139

Identification du gène (NCBI):

16398

Nom complet:

integrin alpha 2

Méthode de purification:

Purification par affinité

CloneNo.:

HMA2

Excitation/Emission maxima
wavelengths:

495 nm / 524 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

souris

Informations générales

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ mouse splenocytes were surface stained with APC Anti-Mouse NK1.1 (CD161) (APC-65138, Clone: PK136) and 0.5 ug FITC Plus Anti-Mouse CD49b (FITC-65245, Clone: HMA 2) or 0.5 ug FITC Plus Armenian Hamster IgG. Cells were not fixed.