

À des fins de recherche uniquement

Anticorps Monoclonal anti-GL7

Numéro de catalogue: **FITC-65261**



Informations de base

Numéro de catalogue:

FITC-65261

Taille:

100ug , 0.5 mg/ml

Hôte:

Rat

Isotype:

IgM, kappa

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

104230

Nom complet:

lymphocyte antigen 77

Méthode de purification:

Purification par affinité

CloneNo.:

GL-7

Excitation/Emission maxima
wavelengths:

495 nm / 524 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

Humain, souris

Informations générales

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ unstimulated or anti-CD3/CD28 stimulated mouse splenocytes were surface co-stained with CoraLite® Plus 647 Anti-Mouse CD3 (17A2) and 0.5 ug FITC Plus Anti-Mouse GL7 (FITC-65261, Clone:GL-7) or 0.5 ug Isotype Control. Cells were not fixed.