

À des fins de recherche uniquement

Anticorps Monoclonal anti-Human IgG Kappa chain



Numéro de catalogue: **HRP-67761**

Informations de base

Numéro de catalogue:

HRP-67761

Taille:

100ul, Concentration: 1000 µg/ml by Nanodrop;

Hôte:

Mouse

Isotype:

IgG1

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

Nom complet:

MW calculé

155 kDa

MW observés:

~25 kDa

Méthode de purification:

Purification par protéine G

CloneNo.:

4D11F7

Dilutions recommandées:

WB 1:1000-1:4000

Applications

Applications testées:

WB

Spécificité de l'espèce:

Humain

Contrôles positifs:

WB : plasma humain,

Informations générales

This antibody detects the Kappa chain of human IgG,

Stockage

Stockage:

Stocker à -20 °C. Éviter toute exposition à la lumière.

Tampon de stockage:

PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



human plasma were subjected to SDS PAGE followed by western blot with HRP-67761 (Human IgG Kappa chain antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.