


# GSK3 SIGNALING PATHWAY

www.ptglab.com



## Related Products

Legend	
—	Association
→	Directly Activates
⊣	Inhibits
Ub	Ubiquitin
P	Phosphorylates

Antibody Name	Catalog Number	Type	Applications
AKT	60203-2-Ig	Mouse Mono	ELISA, IF, IHC, IP, WB
AKT1	10176-2-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
AKT-phospho-S473	66444-1-Ig	Mouse Mono	ELISA, IHC, WB
APC	19782-1-AP	Rabbit Poly	ELISA, IHC
AXIN1	16541-1-AP	Rabbit Poly	ELISA, WB
beta-Catenin	51067-2-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
beta-Catenin	66379-1-Ig	Mouse Mono	ELISA, IHC, WB
Casein Kinase 1 delta	14388-1-AP	Rabbit Poly	ELISA, IF, IP, WB
Casein Kinase 1 epsilon	11230-1-AP	Rabbit Poly	ELISA, IHC, WB
CTBP1	10972-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
eIF2A	11233-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
ERK1/2	11257-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
ERK1/2	66192-1-Ig	Mouse Mono	ELISA, IF, IHC, IP, WB

Antibody Name	Catalog Number	Type	Applications
GBP1	15303-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
GSK3A	15419-1-AP	Rabbit Poly	ELISA, IHC, WB
GSK3B	22104-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
PI3K p110 (alpha)	20583-1-AP	Rabbit Poly	ELISA, FC, IHC, WB
PI3K p110 (beta)	20584-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
PI3K p110 (gamma)	20662-1-AP	Rabbit Poly	ELISA, IP
PP2A	10321-1-AP	Rabbit Poly	ELISA, IHC, WB
PTEN	22054-1-AP	Rabbit Poly	ELISA, IF, IHC, WB

\* This number shows the amount of times our antibody has been cited in a publication.  
 \* More validation images available on our website.

