


# cAMP PATHWAY

www.ptglab.com



## Related Products

Legend	Association
→	Directly Activates
⊣	Inhibits
P	Phosphorylates

Antibody Name	Catalog Number	Type	Applications
14-3-3	14503-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
14-3-3	66061-1-Ig	Mouse Mono	ELISA, FC, IF, IHC, WB
ADCY1	55067-3-AP	Rabbit Poly	ELISA, WB
ADCY3	19492-3-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
AKT1	10176-2-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
AKT1/3	60203-1-Ig	Mouse Mono	ELISA, FC, IF, IHC, IP, WB
AKT-phospho-S473	66444-1-Ig	Mouse Mono	ELISA, IHC, WB
Alpha adducin	10791-1-AP	Rabbit Poly	ELISA, IHC, WB
ATF1	11946-1-AP	Rabbit Poly	ELISA, IF, IP, WB
Beta adducin	14640-1-AP	Rabbit Poly	ELISA, IF, WB
beta Arrestin 1	15361-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
CAMK2	13730-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
CaMKII alpha-Specific	20666-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB

Antibody Name	Catalog Number	Type	Applications
CaMKIV	13263-1-AP	Rabbit Poly	ELISA, IHC, WB
cAMP	12009-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
CREB1	12208-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
CREM	12131-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
DARPP32	10748-1-AP	Rabbit Poly	ELISA, IHC, WB
ERK1/2	11257-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
Gamma adducin	17585-1-AP	Rabbit Poly	ELISA, IP, WB
Glucagon	15954-1-AP	Rabbit Poly	ELISA, IHC, WB
GSK3A	13419-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
GSK3B	22104-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
MEK1-Specific	51080-1-AP	Rabbit Poly	ELISA, IHC
MTOR	20657-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
Myosin 2a	66212-1-Ig	Mouse Mono	ELISA, IHC, WB

Antibody Name	Catalog Number	Type	Applications
NFATC2	22023-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
NF-κB	14220-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
p38	14064-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
p70S6K	14485-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
PVK2	17592-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
RAF1	26863-1-AP	Rabbit Poly	ELISA, IHC, WB
RHOA	10749-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
SRC	11097-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
UCN	21951-1-AP	Rabbit Poly	ELISA, IHC
UCN2	12198-1-AP	Rabbit Poly	ELISA, IHC
Urocortin	19731-1-AP	Rabbit Poly	ELISA, WB

This number shows the amount of times our antibody has been cited in a publication.