



FlexAble Antibody Labeling Kits: Fast & Easy

FlexAble

Any antibody. Any color. Any time.



FlexAble is a novel antibody labeling kit for immunofluorescence, immunohistochemistry, western blot and flow cytometry. It uses an affinity linker to conjugate fluorochromes, enzymes and molecules in any buffer condition.

- **Fast and easy 2-step protocol**
- **Ready to use in 10 minutes**
- **No buffer exchange required**

SPECIFICITIES

Rabbit IgG
Mouse IgG1
Mouse IgG2a
Mouse IgG2b
Human IgG

CONJUGATES

CoraLite® Plus 405
FITC Plus
CoraLite® Plus 488
CoraLite® Plus 555
CoraLite® Plus 647
CoraLite® Plus 750
HRP
Biotin

