



FlexAble Antibody Labeling Kits for Multiplexing

HeLa cells stained with 3 rabbit antibodies: anti-TDP-43 (10782-2-AP) + FlexAble CoraLite® 488 Kit (KFA001, green), anti-TOM20 (11802-1-AP) + FlexAble CoraLite® Plus 555 Kit (KFA002, red), anti-Lamin B1 (12987-1-AP) + FlexAble CoraLite® Plus 647 Kit (cyan, KFA003).

ptglab.com

FlexAble

Any antibody. Any color. Any time.



FlexAble is a novel antibody labeling kit that simplifies multiplexing in immunofluorescence, immunohistochemistry, western blot and flow cytometry. Antibodies are labeled and ready for use in 10 minutes.

- **Simplify multiplex staining with antibodies of the same species**
- **Label antibodies from different species in the same tube**
- **Multiplex your FlexAble-labeled antibodies with:**
 - Other antibodies labeled with a FlexAble Kit
 - Pre-conjugated antibodies from any supplier
 - Primary and secondary antibodies

SPECIFICITIES

Rabbit IgG
Mouse IgG1
Mouse IgG2a
Mouse IgG2b
Human IgG

CONJUGATES

CoraLite® Plus 405
FITC Plus
CoraLite® Plus 488
CoraLite® Plus 555
CoraLite® Plus 647
CoraLite® Plus 750
HRP
Biotin

