

DYKDDDDK Fab-Trap™


A premium resin for high purity immunoprecipitation of Flag®-tagged proteins.

- Lowest background in class
- Less contaminating antibody chains
- Recombinantly expressed Fab-fragment derived from M2 clone
- Compatible with harsh buffer conditions
- Ready-to-use

For more information & free sample, visit our website.



DYKDDDDK Fab-Trap™ consisting of an anti-Flag® Fab-fragment (blue) conjugated to Agarose beads (grey).



Immunoprecipitation of Flag® fusion protein with DYKDDDDK Fab-Trap™, elution with 3xDYKDDDDK-peptide. I: Input, FT: Flow-Through, E: Elution.

DYKDDDDK reagents for Immunoprecipitation & Western blot

- DYKDDDDK Fab-Trap™ Agarose
- 3xDYKDDDDK-peptide
- DYKDDDDK tag Polyclonal antibody
- DYKDDDDK Immunoprecipitation Starter Pack

For more information & free sample visit our
website or contact:

ChromoTek
E: order@chromotek.com

Proteintech Group
E: proteintech@ptglab.com



SCAN ME