

For Research Use Only

Fas Ligand Polyclonal antibody

Catalog Number: 11322-1-AP **2 Publications**



Basic Information

Catalog Number: 11322-1-AP	GenBank Accession Number: BC017502	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 356	
Source: Rabbit	UNIPROT ID: P48023	
Isotype: IgG	Full Name: Fas ligand (TNF superfamily, member 6)	
Immunogen Catalog Number: AG1816	Calculated MW: 281 aa, 31 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human

Cited Species:
human, pig

Notable Publications

Author	Pubmed ID	Journal	Application
Ben-Cheng Lin	31954379	J Biochem Mol Toxicol	WB
Jian Li	38508351	Comp Biochem Physiol C Toxicol Pharmacol	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data