

For Research Use Only

ATCAY Polyclonal antibody

Catalog Number: 13257-2-AP



Basic Information

Catalog Number:

13257-2-AP

Size:

150ul , Concentration: 280 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3658

GenBank Accession Number:

BC026217

GeneID (NCBI):

85300

UNIPROT ID:

Q86WG3

Full Name:

ataxia, cerebellar, Cayman type

Calculated MW:

371 aa, 42 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 13257-2-AP (ATCAY antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.