

For Research Use Only

# TAB1 Polyclonal antibody

Catalog Number: 14819-1-AP

Featured Product

5 Publications



## Basic Information

<b>Catalog Number:</b> 14819-1-AP	<b>GenBank Accession Number:</b> BC050554	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul, Concentration: 380 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 10454	
<b>Source:</b> Rabbit	<b>Full Name:</b> mitogen-activated protein kinase kinase kinase 7 interacting protein 1	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 55 kDa	
<b>Immunogen Catalog Number:</b> AG6590		

## Applications

**Tested Applications:**  
ELISA

**Cited Applications:**  
CoIP, WB

**Species Specificity:**  
human, mouse, rat

**Cited Species:**  
human, sheep

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Zhongjia Jiang	28553017	Mediators Inflamm	WB
Lin Yuan	35352799	Cardiovasc Res	WB, CoIP
Kun-Lin Chen	31285783	Oxid Med Cell Longev	WB

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data