

For Research Use Only

UBE2Q1 Polyclonal antibody

Catalog Number: 16072-1-AP



Basic Information

Catalog Number:

16072-1-AP

Size:

150ul , Concentration: 113 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9028

GenBank Accession Number:

BC009286

GeneID (NCBI):

55585

UNIPROT ID:

Q7Z7E8

Full Name:

ubiquitin-conjugating enzyme E2Q family member 1

Calculated MW:

422 aa, 46 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-UBE2Q1 (IP:16072-1-AP, 4ug;
Detection:16072-1-AP 1:600) with HEK-293 cells
lysate 1440 ug.