

For Research Use Only

KATNAL2 Polyclonal antibody, PBS Only

Catalog Number: 17074-1-PBS



Basic Information

Catalog Number:

17074-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10950

GenBank Accession Number:

BC034999

GeneID (NCBI):

83473

UNIPROT ID:

Q8IYT4

Full Name:

katanin p60 subunit A-like 2

Calculated MW:

466 aa, 53 kDa

Observed MW:

40-80 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IF/ICC, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using KATNAL2 antibody (17074-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 17074-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed NIH/3T3 cells using KATNAL2 antibody (17074-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 17074-1-PBS in a different storage buffer formulation.