

For Research Use Only

# TLX3 Polyclonal antibody

Catalog Number: 17454-1-AP



## Basic Information

**Catalog Number:**

17454-1-AP

**Size:**

150ul , Concentration: 500 µg/ml by Nanodrop and 300 µg/ml by Bradford method using BSA as the standard;

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG10612

**GenBank Accession Number:**

BC017291

**GeneID (NCBI):**

30012

**UNIPROT ID:**

O43711

**Full Name:**

T-cell leukemia homeobox 3

**Calculated MW:**

291 aa, 32 kDa

**Observed MW:**

32 kDa

**Purification Method:**

Antigen affinity purification

**Recommended Dilutions:**

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

## Applications

**Tested Applications:**

IP, ELISA

**Species Specificity:**

human, mouse, rat

**Positive Controls:**

IP : Jurkat cells,

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



IP result of anti-TLX3 (IP:17454-1-AP, 4ug;  
Detection:17454-1-AP 1:300) with Jurkat cells  
lysate 2800ug.