

For Research Use Only

Integrin alpha 10 Specific Polyclonal antibody



Catalog Number: 19737-1-AP

Basic Information

Catalog Number:

19737-1-AP

Size:

150ul , Concentration: 287 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_003637

GeneID (NCBI):

8515

UNIPROT ID:

O75578

Full Name:

integrin, alpha 10

Calculated MW:

128 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

ITGA10, also named as UNQ468 and PRO827, belongs to the integrin alpha chain family. Integrin alpha-10/beta-1 is a receptor for collagen. This antibody is specific to ITGA10.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data