

For Research Use Only

FOXJ3 Polyclonal antibody, PBS Only

Catalog Number:19751-1-PBS



Basic Information

Catalog Number:

19751-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_014947

GeneID (NCBI):

22887

UNIPROT ID:

Q9UPW0

Full Name:

forkhead box J3

Calculated MW:

69 kDa

Observed MW:

69 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

FOXJ3, also named as KIAA1041, contains 1 fork-head DNA-binding domain. FOXJ3 is a transcriptional activator.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human liver tissue were subjected to SDS PAGE followed by western blot with 19751-1-AP (FOXJ3 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 19751-1-PBS in a different storage buffer formulation.