

For Research Use Only

ARMCX1 Polyclonal antibody, PBS Only

Catalog Number: 20193-1-PBS

Featured Product



Basic Information

Catalog Number:

20193-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14120

GenBank Accession Number:

BC002691

GeneID (NCBI):

51309

UNIPROT ID:

Q9P291

Full Name:

armadillo repeat containing, X-linked 1

Calculated MW:

453 aa, 49 kDa

Observed MW:

49-55 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse ovary tissue were subjected to SDS PAGE followed by western blot with 20193-1-AP (ARMCX1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20193-1-PBS in a different storage buffer formulation.