

For Research Use Only

POTEC-Specific Polyclonal antibody

Catalog Number: 20204-1-AP



Basic Information

Catalog Number: 20204-1-AP	GenBank Accession Number: NM_001137671	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 153 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 388468	
Source: Rabbit	UNIPROT ID: B2RU33	
Isotype: IgG	Full Name: POTE ankyrin domain family, member C	
	Calculated MW: 61 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

POTEC, also named as A26B2 and POTE18, belongs to the POTE family. POTE has an important signaling function in the reproductive system. This antibody is specific to POTEC.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data