

For Research Use Only

METTL20 Polyclonal antibody

Catalog Number: 24448-1-AP



Basic Information

Catalog Number: 24448-1-AP	GenBank Accession Number: BC039107	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 900 µg/ml by Nanodrop and 540 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 254013	
Source: Rabbit	UNIPROT ID: Q8IXQ9	
Isotype: IgG	Full Name: chromosome 12 open reading frame 72	
Immunogen Catalog Number: AG19798	Calculated MW: 262 aa, 29 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data