

For Research Use Only

FAM47B Polyclonal antibody

Catalog Number: 24467-1-AP



Basic Information

Catalog Number: 24467-1-AP	GenBank Accession Number: BC035026	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 750 ug/ml by Nanodrop;	GeneID (NCBI): 170062	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q8NA70	
Isotype: IgG	Full Name: family with sequence similarity 47, member B	
Immunogen Catalog Number: AG21562	Calculated MW: 645 aa, 74 kDa	
	Observed MW: 70 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse testis tissue, rat testis tissue
Species Specificity: human, mouse, rat	

Background Information

FAM47B(Family with sequence similarity 47 member B) belongs to the FAM47 family and is expressed in testis.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24467-1-AP (FAM47B antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.