

For Research Use Only

TRIM40 Polyclonal antibody

Catalog Number: 24526-1-AP



Basic Information

Catalog Number: 24526-1-AP	GenBank Accession Number: BC060785	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 450 ug/ml by Nanodrop;	GeneID (NCBI): 135644	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	UNIPROT ID: Q6P9F5	
Isotype: IgG	Full Name: tripartite motif-containing 40	
Immunogen Catalog Number: AG21649	Calculated MW: 258 aa, 29 kDa	
	Observed MW: 29-35 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse liver tissue, rat liver tissue
Species Specificity: mouse, rat	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24526-1-AP (TRIM40 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.