

For Research Use Only

LRRC4 Polyclonal antibody

Catalog Number: 24834-1-AP



Basic Information

Catalog Number: 24834-1-AP	GenBank Accession Number: BC111561	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 260 ug/ml by Nanodrop;	GeneID (NCBI): 64101	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	UNIPROT ID: Q9HBW1	
Isotype: IgG	Full Name: leucine rich repeat containing 4	
Immunogen Catalog Number: AG20376	Calculated MW: 653 aa, 73 kDa	
	Observed MW: 73 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse liver tissue,
Species Specificity: human, mouse	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24834-1-AP (LRRC4 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.