

For Research Use Only

PLEKHF2 Polyclonal antibody

Catalog Number: 25424-1-AP



Basic Information

Catalog Number:

25424-1-AP

Size:

150ul , Concentration: 450 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22285

GenBank Accession Number:

BC011806

GeneID (NCBI):

79666

UNIPROT ID:

Q9H8W4

Full Name:

pleckstrin homology domain containing, family F (with FYVE domain) member 2

Calculated MW:

249 aa, 28 kDa

Observed MW:

22-27 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, rat brain tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25424-1-AP (PLEKHF2 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.