

For Research Use Only

TTC38 Polyclonal antibody

Catalog Number: 25452-1-AP



Basic Information

Catalog Number:

25452-1-AP

Size:

150ul , Concentration: 450 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22070

GenBank Accession Number:

BC018918

GeneID (NCBI):

55020

UNIPROT ID:

Q5R3I4

Full Name:

tetratricopeptide repeat domain 38

Calculated MW:

469 aa, 53 kDa

Observed MW:

48-53 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:5000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HepG2 cells, L02 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25452-1-AP (TTC38 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.