

For Research Use Only

MORG1 Polyclonal antibody

Catalog Number: 27244-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number: 27244-1-AP	GenBank Accession Number: BC005870	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 500 µg/ml by Nanodrop and 300 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 84292	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: mitogen-activated protein kinase organizer 1	
Isotype: IgG	Calculated MW: 315 aa, 34 kDa	
Immunogen Catalog Number: AG15141	Observed MW: 30-35 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse testis tissue,
Cited Applications: IF, WB	
Species Specificity: human, mouse	
Cited Species: mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Guoyao Ou	37210913	Int Immunopharmacol	WB,IF

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 27244-1-AP (MORG1 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.