

For Research Use Only

Phospho-SMAD2 (Ser465/467) Polyclonal antibody



Catalog Number: 29129-1-AP **1 Publications**

Basic Information

Catalog Number:

29129-1-AP

Size:

100ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC014840

GeneID (NCBI):

4087

UNIPROT ID:

Q15796

Full Name:

SMAD family member 2

Calculated MW:

467 aa, 52 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

Human

Cited Species:

human, mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Jinjin Jiang	35512472	Ecotoxicol Environ Saf	WB

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data