

For Research Use Only

ZHX3 Polyclonal antibody

Catalog Number: 29397-1-AP **1 Publications**



Basic Information

Catalog Number: 29397-1-AP	GenBank Accession Number: NM_015035	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 1100 ug/ml by Nanodrop;	GeneID (NCBI): 23051	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	UNIPROT ID: Q9H4I2	
Isotype: IgG	Full Name: zinc fingers and homeoboxes 3	
Immunogen Catalog Number: AG30023	Calculated MW: 105 kDa	
	Observed MW: 130 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HEK-293 cells, Jurkat cells, MCF-7 cells, HeLa cells
Cited Applications: WB, IF	
Species Specificity: Human	
Cited Species: human	

Background Information

ZHX3 is a member of the zinc fingers and homeoboxes (ZHX) gene family. ZHX3 is a ubiquitous transcriptional repressor expressed in human cells and tissues.

Notable Publications

Author	Pubmed ID	Journal	Application
Dongwei Yao	39267687	Am J Cancer Res	WB,IF

Storage


Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29397-1-AP (ZHX3 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.