

For Research Use Only

# STIL Polyclonal antibody

Catalog Number: 30504-1-AP



## Basic Information

Catalog Number:

30504-1-AP

Size:

150ul , Concentration: 1400 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG32611

GenBank Accession Number:

BC126223

GeneID (NCBI):

6491

UNIPROT ID:

Q15468

Full Name:

SCL/TAL1 interrupting locus

Calculated MW:

1287 aa, 143 kDa

Observed MW:

150 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : Jurkat cells, K-562 cells, THP-1 cells

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30504-1-AP (STIL antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.