

For Research Use Only

ACIN1-L specific Polyclonal antibody

Catalog Number:30522-1-AP



Basic Information

Catalog Number: 30522-1-AP	GenBank Accession Number: BC140805	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 350 ug/ml by Nanodrop;	GeneID (NCBI): 22985	Recommended Dilutions: WB 1:500-1:1000 IF-P 1:50-1:500
Source: Rabbit	UNIPROT ID: Q9UKV3	
Isotype: IgG	Full Name: apoptotic chromatin condensation inducer 1	
Immunogen Catalog Number: AG30506	Calculated MW: 1341 aa, 152 kDa	
	Observed MW: 240 kDa	

Applications

Tested Applications: WB, IF-P, ELISA	Positive Controls: WB : HEK-293 cells, HeLa cells, HepG2 cells
Species Specificity: Human, mouse	IF-P : mouse cerebellum tissue,

Background Information

Apoptosis is defined by several morphologic nuclear changes, including chromatin condensation and nuclear fragmentation. ACIN1 is a nuclear protein that induces apoptotic chromatin condensation after activation by caspase-3, without inducing DNA fragmentation. This protein has also been shown to be a component of a splicing-dependent multiprotein exon junction complex (EJC) that is deposited at splice junctions on mRNAs, as a consequence of pre-mRNA splicing. There are three isoforms of human acinus (the L, S and S' forms) shows at about 80-85 kDa and 220-250 kDa by Western Blot, which were probably generated by alternative splicing (PMID: 10490026). 30522-1-AP specifically targets at the ACIN1-L, while another antibody 23937-1-AP can detect all three isoforms of ACIN1.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30522-1-AP (ACIN1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed mouse cerebellum tissue using ACIN1-L specific antibody (30522-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).