

For Research Use Only

NPNT Polyclonal antibody

Catalog Number: 30836-1-AP



Basic Information

Catalog Number:

30836-1-AP

Size:

150ul , Concentration: 500 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34341

GenBank Accession Number:

BC166654

GeneID (NCBI):

255743

UNIPROT ID:

Q6UXI9

Full Name:

nephrinectin

Calculated MW:

62 kDa

Observed MW:

62-70kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse , rat

Positive Controls:

WB : HEK-293 cells, HepG2 cells, mouse kidney tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30836-1-AP (NPNT antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.