

For Research Use Only

TYK2 Polyclonal antibody

Catalog Number: 30838-1-AP



Basic Information

Catalog Number:

30838-1-AP

Size:

150ul , Concentration: 400 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28457

GenBank Accession Number:

BC014243

GeneID (NCBI):

7297

UNIPROT ID:

P29597

Full Name:

tyrosine kinase 2

Calculated MW:

1187 aa, 134 kDa

Observed MW:

134 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, HeLa cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30838-1-AP (TYK2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.