

For Research Use Only

HDGFL1 Polyclonal antibody

Catalog Number: 31923-1-AP



Basic Information

Catalog Number: 31923-1-AP	GenBank Accession Number: NM_138574	Purification Method: Antigen affinity Purification
Size: 150ul , Concentration: 550 ug/ml by Nanodrop;	GeneID (NCBI): 154150	Recommended Dilutions: WB 1:1000-1:6000
Source: Rabbit	UNIPROT ID: Q5TGJ6	
Isotype: IgG	Full Name: hepatoma derived growth factor-like 1	
Immunogen Catalog Number: AG36944	Calculated MW: 27 kDa	
	Observed MW: 40 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : human testis tissue,
Species Specificity: human	

Background Information

Hepatoma-derived growth factor-like 1 (HDGFL1, also known as HRP-1, PWWP1, and Hdgfrp1) is a member of the hepatoma-derived growth factor (HDGF) family related to tumorigenesis and the development of cancer (PMID: 18823380; 20509080). It is specifically expressed in the testis and activated in the nucleus (PMID: 22217008).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31923-1-AP (HDGFL1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.