

For Research Use Only

ABCC2 Polyclonal antibody

Catalog Number: 59012-1-Ig



Basic Information

Catalog Number:

59012-1-Ig

Size:

150ul , Concentration: 1000 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC136419

GeneID (NCBI):

1244

UNIPROT ID:

Q92887

Full Name:

ATP-binding cassette, sub-family C (CFTR/MRP), member 2

Calculated MW:

1545 aa, 174 kDa

Purification Method:

Protein G purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data