

For Research Use Only

BAX Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60267-2-PBS



Basic Information

Catalog Number: 60267-2-PBS	GenBank Accession Number: BC014175	Purification Method: Protein A Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 581	CloneNo.: 5G3G7
Source: Mouse	UNIPROT ID: Q07812	
Isotype: IgG2a	Full Name: BCL2-associated X protein	
Immunogen Catalog Number: AG21068	Calculated MW: 21 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

60267-2-PBS targets BAX as part of a matched antibody pair:

MP50504-1: 60267-2-PBS capture and 60267-3-PBS detection (validated in Cytometric bead array)

MP50504-2: 60267-4-PBS capture and 60267-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50504-1, BAX Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60267-2-PBS. Detection antibody: 60267-3-PBS. Standard:Ag21068. Range: 0.391-25 ng/mL



Cytometric bead array standard curve of MP50504-2, BAX Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60267-4-PBS. Detection antibody: 60267-2-PBS. Standard:Ag21068. Range: 0.098-100 ng/mL