

For Research Use Only

FGFR2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60403-1-PBS



Basic Information

Catalog Number: 60403-1-PBS	GenBank Accession Number: 2263	Purification Method: Protein A Magarose purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 2263	CloneNo.: 3G11A2
Source: Mouse	UNIPROT ID: P21802-3	
Isotype: IgG2a	Full Name: fibroblast growth factor receptor 2	
	Calculated MW: 92kd	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

60403-1-PBS targets FGFR2 as part of a matched antibody pair.

MP50495-1: 60403-1-PBS capture and 60403-2-PBS detection (validated in Cytometric bead array)

MP50495-3: 60403-1-PBS capture and 60403-5-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50495-1, FGFR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60403-1-PBS. Detection antibody: 60403-2-PBS. Standard:Eg0175. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50495-3, FGFR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60403-1-PBS. Detection antibody: 60403-5-PBS. Standard:Eg0175. Range: 0.391-25 ng/mL