

For Research Use Only

EVI1 Monoclonal antibody, PBS Only

Catalog Number:60413-1-PBS



Basic Information

Catalog Number: 60413-1-PBS	GenBank Accession Number: NM_001105077	Purification Method: Protein G purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 2122	CloneNo.: 1A12A6
Source: Mouse	UNIPROT ID: Q03112	
Isotype: IgG1	Full Name: ecotropic viral integration site 1	
Immunogen Catalog Number: AG34143	Calculated MW: 138 kDa	
	Observed MW: 140-150 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MKN-45 cells were subjected to SDS PAGE followed by western blot with 60413-1-Ig (EVI1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 60413-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 60413-1-Ig (EVI1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 60413-1-PBS in a different storage buffer formulation.