

For Research Use Only

# CLTA Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60427-3-PBS



## Basic Information

<b>Catalog Number:</b> 60427-3-PBS	<b>GenBank Accession Number:</b> BC019287	<b>Purification Method:</b> Protein A Magarose purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 1211	<b>CloneNo.:</b> 2C5F10
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P09496	
<b>Isotype:</b> IgG2a	<b>Full Name:</b> clathrin, light chain (Lca)	
<b>Immunogen Catalog Number:</b> AG1299	<b>Calculated MW:</b> 27 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

60427-3-PBS targets CLTA as part of a matched antibody pair:

MP50560-4: 60427-3-PBS capture and 60427-4-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50560-4, CLTA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60427-3-PBS. Detection antibody: 60427-4-PBS. Standard: Ag1299. Range: 1.563-100 ng/mL.