

For Research Use Only

# CHCHD10 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60475-2-PBS



## Basic Information

**Catalog Number:**

60475-2-PBS

**Size:**

100ug, Concentration: 1 mg/ml by Nanodrop;

**Source:**

Mouse

**Isotype:**

IgG1

**Immunogen Catalog Number:**

AG22598

**GenBank Accession Number:**

BC065232

**GeneID (NCBI):**

400916

**UNIPROT ID:**

Q8WYQ3

**Full Name:**

coiled-coil-helix-coiled-coil-helix domain containing 10

**Calculated MW:**

14 kDa

**Purification Method:**

Protein G Magarose purification

**CloneNo.:**

1B5E2

## Applications

**Tested Applications:**

Cytometric bead array, Indirect ELISA

**Species Specificity:**

human

## Product Information

60475-2-PBS targets CHCHD10 as part of a matched antibody pair:

MP50651-1: 60475-1-PBS capture and 60475-2-PBS detection (validated in Cytometric bead array)

MP50651-4: 60475-4-PBS capture and 60475-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50651-1, CHCHD10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60475-1-PBS. Detection antibody: 60475-2-PBS. Standard: Ag22598. Range: 0.098-25 ng/mL



Cytometric bead array standard curve of MP50651-4, CHCHD10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60475-4-PBS. Detection antibody: 60475-2-PBS. Standard: Ag22598. Range: 1.563-100 ng/mL