

For Research Use Only

FSTL3/FLRG Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60493-3-PBS



Basic Information

Catalog Number:

60493-3-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM_005860.3

GeneID (NCBI):

10272

UNIPROT ID:

O95633-1

Full Name:

follicle-stimulating-like 3 (secreted glycoprotein)

Calculated MW:

28 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3G10B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Product Information

60493-3-PBS targets FSTL3/FLRG as part of a matched antibody pair:

MP50689-2: 60493-3-PBS capture and 60493-4-PBS detection (validated in Cytometric bead array)

MP50689-3: 60493-3-PBS capture and 60493-5-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50689-2, FSTL3/FLRG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60493-3-PBS. Detection antibody: 60493-4-PBS. Standard:Eg0949. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50689-3, FSTL3/FLRG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60493-3-PBS. Detection antibody: 60493-5-PBS. Standard:null. Range: null.