

For Research Use Only

Ubiquilin 2 Monoclonal antibody, PBS Only

Catalog Number: 60495-1-PBS



Basic Information

Catalog Number:

60495-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG20251

GenBank Accession Number:

BC069237

GeneID (NCBI):

29978

UNIPROT ID:

Q9UHD9

Full Name:

ubiquilin 2

Calculated MW:

624 aa, 66 kDa

Observed MW:

66 kDa

Purification Method:

Protein G purification

CloneNo.:

1A6D11

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human, rat, pig, rabbit

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 60495-1-Ig (Ubiquilin 2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 60495-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed A431 cells using Ubiquilin 2 antibody (60495-1-Ig, Clone: 1A6D11) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Mouse IgG(H+L) (SA00013-1). This data was developed using the same antibody clone with 60495-1-PBS in a different storage buffer formulation.