

For Research Use Only

SMCR7,MID49 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60512-3-PBS



Basic Information

Catalog Number:

60512-3-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG9503

GenBank Accession Number:

BC014973

GeneID (NCBI):

125170

UNIPROT ID:

Q96C03

Full Name:

Smith-Magenis syndrome
chromosome region, candidate 7

Calculated MW:

454 aa, 49 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1F8G8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Product Information

60512-3-PBS targets SMCR7,MID49 as part of a matched antibody pair:

MP50716-2: 60512-3-PBS capture and 60512-4-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50716-2, SMCR7,MID49 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60512-3-PBS. Detection antibody: 60512-4-PBS. Standard:Ag9503. Range: 0.391-25 ng/mL.