

For Research Use Only

# HSD17B2 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60587-1-PBS



## Basic Information

<b>Catalog Number:</b> 60587-1-PBS	<b>GenBank Accession Number:</b> BC009581	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 100ug, Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 3294	<b>CloneNo.:</b> 2C1A2
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P37059	
<b>Isotype:</b> IgG1	<b>Full Name:</b> hydroxysteroid (17-beta) dehydrogenase 2	
<b>Immunogen Catalog Number:</b> AG1411	<b>Calculated MW:</b> 43 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

60587-1-PBS targets HSD17B2 as part of a matched antibody pair:

MP50836-1: 60587-1-PBS capture and 60587-2-PBS detection (validated in Cytometric bead array)

MP50836-2: 60587-1-PBS capture and 60587-3-PBS detection (validated in Cytometric bead array)

MP50836-3: 60587-1-PBS capture and 60587-4-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data



Cytometric bead array standard curve of MP50836-1, HSD17B2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60587-1-PBS. Detection antibody: 60587-2-PBS. Standard:Ag1411. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP50836-2, HSD17B2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60587-1-PBS. Detection antibody: 60587-3-PBS. Standard:Ag1411. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50836-3, HSD17B2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60587-1-PBS. Detection antibody: 60587-4-PBS. Standard:Ag1411. Range: 0.781-100 ng/mL