

For Research Use Only

SERCA2/ATP2A2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60593-2-PBS



Basic Information

Catalog Number: 60593-2-PBS	GenBank Accession Number: BC035588	Purification Method: Protein G Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 488	CloneNo.: 1F3G1
Source: Mouse	UNIPROT ID: P16615	
Isotype: IgG1	Full Name: ATPase, Ca ⁺⁺ transporting, cardiac muscle, slow twitch 2	
Immunogen Catalog Number: AG26307	Calculated MW: 1042 aa, 115 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

60593-2-PBS targets SERCA2/ATP2A2 as part of a matched antibody pair.

MP50843-1: 60593-1-PBS capture and 60593-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50843-1, SERCA2,ATP2A2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60593-1-PBS. Detection antibody: 60593-2-PBS. Standard:Ag5052. Range: 0.098-25 ng/mL.