

For Research Use Only

# TLR3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60598-1-PBS



## Basic Information

<b>Catalog Number:</b> 60598-1-PBS	<b>GenBank Accession Number:</b> BC094737	<b>Purification Method:</b> Protein G Magarose purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 7098	<b>CloneNo.:</b> 2A1A7
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> O15455	
<b>Isotype:</b> IgG1	<b>Full Name:</b> toll-like receptor 3	
<b>Immunogen Catalog Number:</b> AG12163	<b>Calculated MW:</b> 904 aa, 104 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

60598-1-PBS targets TLR3 as part of a matched antibody pair:

MP50854-1: 60598-1-PBS capture and 60598-2-PBS detection (validated in Cytometric bead array)

MP50854-2: 60598-1-PBS capture and 60598-3-PBS detection (validated in Cytometric bead array)

MP50854-3: 60598-1-PBS capture and 60598-4-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data



Cytometric bead array standard curve of MP50854-1, TLR3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60598-1-PBS. Detection antibody: 60598-2-PBS. Standard:Ag12163. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50854-2, TLR3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60598-1-PBS. Detection antibody: 60598-3-PBS. Standard:Ag12163. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50854-3, TLR3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60598-1-PBS. Detection antibody: 60598-4-PBS. Standard:Ag12163. Range: 0.781-50 ng/mL.