

For Research Use Only

# SURVIVIN-Specific Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60760-1-PBS



## Basic Information

Catalog Number:

60760-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG3

GenBank Accession Number:

NM\_001168

GeneID (NCBI):

332

UNIPROT ID:

O15392

Full Name:

baculoviral IAP repeat-containing 5

Calculated MW:

18 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

5G1B2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Product Information

60760-1-PBS targets SURVIVIN-Specific as part of a matched antibody pair:

MP51208-3: 66495-2-PBS capture and 60760-1-PBS detection (validated in Cytometric bead array)

MP51208-4: 66495-4-PBS capture and 60760-1-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51208-3, SURVIVIN Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66495-2-PBS. Detection antibody: 60760-1-PBS. Standard:Ag20958. Range: 3.125-100 ng/mL



Cytometric bead array standard curve of MP51208-4, SURVIVIN Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66495-4-PBS. Detection antibody: 60760-1-PBS. Standard:Ag20958. Range: 3.125-100 ng/mL