

For Research Use Only

CD268/BAFF-R Monoclonal antibody, PBS Only (Detector)

Catalog Number:60763-4-PBS



Basic Information

Catalog Number:

60763-4-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgM

Immunogen Catalog Number:

AG18512

GenBank Accession Number:

BC105123

GeneID (NCBI):

115650

UNIPROT ID:

Q96RJ3

Full Name:

tumor necrosis factor receptor
superfamily, member 13C

Calculated MW:

184 aa, 19 kDa

Purification Method:

Caprylic acid/ammonium sulfate
precipitation

CloneNo.:

2E8E12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Product Information

60763-4-PBS targets CD268/BAFF-R as part of a matched antibody pair.

MP51091-3: 60763-2-PBS capture and 60763-4-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Background Information

B-cell-activating factor receptor (BAFF-R, also known as CD268 and TNFRSF13C) is a member of the tumor necrosis factor receptor (TNFR) superfamily. It is expressed in most malignant B cells, including those in primary central nervous system lymphoma (PCNSL), and is involved in the survival and proliferation of these cells (PMID: 32528875). The interaction of BAFF with BAFF-R promotes NF- κ B activation, which plays a crucial role in B-cell maturation, survival, and costimulation of T-cell activation and proliferation (PMID: 30349534). Additionally, BAFF-R is associated with chronic kidney disease (CKD) by activating NF- κ B signaling in renal tubular epithelial cells (PMID: 38791453).

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51091-3, TNFRSF13C Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60763-2-PBS. Detection antibody: 60763-4-PBS. Standard:Ag18512. Range: 0.391-100 ng/mL