

For Research Use Only

APOA5 Monoclonal antibody

Catalog Number: 65015-1-Ig



Basic Information

Catalog Number: 65015-1-Ig	GenBank Accession Number: BC101787	Purification Method: Antigen affinity purification
Size: 150ul , 500 µg/ml	GeneID (NCBI): 116519	Recommended Dilutions: WB 1:500-1:2000
Source: Mouse	Full Name: apolipoprotein A-V	
Isotype: IgG1	Calculated MW: 366 aa, 41 kDa	
	Observed MW: 39 kDa-41 kDa	

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : Human serum,
Species Specificity: human	

Background Information

Storage

Storage:
Store at 2-8°C.
Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data