

For Research Use Only

# Anti-Mouse CD24 (M1/69)

Catalog Number: 65121-1-Ig



## Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
65121-1-Ig	BC075622	Affinity purification
Size:	GeneID (NCBI):	CloneNo.:
100ug, 500 µg/ml	12484	M1/69
Source:	Full Name:	
Rat	CD24a antigen	
Isotype:		
IgG2b, kappa		

## Applications

Tested Applications:  
FC

Species Specificity:  
Mouse

## Background Information

### Storage


Storage:  
Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.1% sodium azide, pH 7.3.

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> BALB/c mouse splenocytes were surface stained with 0.5 ug Anti-Mouse CD24 (65121-1-Ig, Clone: M1/69) and FITC anti-rat IgG2b Antibody at dilution 1:200 (green), or stained with 0.5 ug isotype control and FITC anti-rat IgG2b Antibody at dilution 1:200 (black). Cells were not fixed.